

# Train and Go Operating Procedure

Garrett Bradshaw, Slade Hicks, Brandon Waldrup, and Kyler Smith

## Purpose

Train and go is a training device designed to help users master the skills required to maneuver a wheelchair. This Device uses virtual reality simulation coupled with a separate operation detection system to help the user achieve this goal. Below is a detailed guide of how to install/operate Train and Go.

## Key Terms/Definitions

*Virtual Reality:* A virtual environment uploaded to the Meta Quest Pro headset.

*Ultrasonic Sensors:* Sensors to be placed along the base of the chair to detect obstacles.

*Microcontroller:* An integrated circuit device that governs a specific operation in an embedded system.


*Printed Circuit Board (PCB):* Circuit board with all the necessary connections for the microcontroller and ultrasonic sensors to operate.

*Inertial Measurement Unit (IMU):* A device that can measure and report specific gravity and angular rate of an object to which it is attached.

*Rumble Motor:* Mechanical device that generates vibrations.

## Installation

- 1) Using the Velcro straps, secure the housing box containing the microcontroller, PCB, and battery to the back of the chair.
  - a. This location should allow easy charging access for the battery.



- 2) Using the Velcro Straps, secure the Meta Quest Pro controller mount to the back of the head rest of the chair.
- 3) Using the provided nuts, secure the ultrasonic sensor mounts along the rail on the sides of the chair.
  - a. There are 6 sensors and 6 sensor mounts. Each mount is designed to face a specific direction (Forward, Adjacent, Backward) on each side of the chair.



- 4) Place the ultrasonic sensors in each sensor mount.
- 5) Connect each wire from the housing box to the corresponding sensor.
  - a. Wires can be secured to the siding of the chair.

- 6) Secure rumble motor inside the back cushion.
  - a. This should be secured in the upper portion of the chair to ensure the user will feel the intended vibration.
  - b. Wiring can be secured to the chair once the rumble motor is set in place.
- 7) Make sure all connections are secured and flip the switch on the housing box.
  - a. Test that each sensor is properly triggering the rumble motor to vibrate at the specified distances in the table below.

Sensor	Detection Range Start (mm)	Detection Range End (mm)
Forward-Left	10	20
Adjacent-Left	10	20
Backward-Left	10	20
Forward-Right	10	20
Adjacent-Right	10	20
Backward-Right	10	20

## Virtual Reality Installation/Setup

Sideload the file TrainandGo.apk onto a VR headset. It is designed for use on Meta Quest Pro. Sideload instructions for a Meta Quest device are as follows:

- 1) Install SideQuest or any sideloading software to a source computer.
  - a. This computer does not need to be powerful enough for VR: it will just transfer files.
  - b. See <https://sidequestvr.com/setup-howto> for installation instructions. Use the “Advanced Installer.”
- 2) Enable developer mode on Meta Quest device.  
<https://help.arborxr.com/en/articles/6333136-developer-mode-usb-debugging-on-meta-quest-devices>
- 3) Connect Meta Quest Pro to the source computer.
- 4) Download TrainandGo.apk to the source computer.
- 5) Load TrainandGo.apk onto the Meta Quest device using SideQuest.
  - a. Drag TrainandGo.apk onto the green bubble in the SideQuest window.
- 6) When SideQuest informs you that the transfer is complete, disconnect the headset.

The remaining VR Instructions are as follows:

- 7) Using the VR Headset, map the boundaries of the room.
- 8) Place the right-hand VR controller in the controller mount that was installed on the back of the wheelchair.



- 9) Use the left-hand VR controller to navigate to the Unknown Sources section of the Apps menu.
- 10) Open TrainandGo.apk.
- 11) Press the trigger on the previously mounted right-hand controller to ensure the controller is “awake” and track the wheelchair movement properly.
- 12) (Optional) Screen share the headset display to a screen for viewing.
  - a. This can be done through SideQuest or Meta’s screen share service. See Links Below: